

## **CUTBACK BITUMEN MC-3000 SPECIFICATIONS**

Parameters	Min	Max	Test Method	
Kinematic viscosity at 60°C mm^2/s	3000	6000	ASTM D2170	
Flash point (tag open-cup),°c	66	-	ASTM D3143	
Distillate test: distillate, volume percent of total distillate at 360°C:				
To 225°C	- (4	-	ASTM D402	
То 260°С	-	15	ASTM D402	
To 316°C	15	75	ASTM D402	
Residue from distillation to 360°c, percent volume by difference	80	-	ASTM D402	
Test on residue from distillation:	est on residue from distillation:			
Penetration at 25°c	120	250	ASTM D5	
Ductility at 25°c	100	-BI	ASTM D113	
Solubility in trichloroethylene, %	99	-	ASTM 2042	
Water, %	-	0.2	ASTM D95	