

CUTBACK BITUMEN SC-250 SPECIFICATIONS

Parameters	Min	Max	Test Method
Kinematic viscosity at 60°C mm²/s	250	500	ASTM D2170
Flash point (tag open-cup), °c	79	-	ASTM D3143
Distillate test: distillate, volume percent of total distillate at 360°C:			
Total distillation up to 360°C, volume%	4	20	ASTM D402
Solubility in Trichloroethylene, %	99	-	ASTM D402
Viscosity at 60°C	800	10000	ASTM D2171
Test on residue from distillation:			
Penetration at 25°C	60	-	ASTM D5
Ductility at 25°C	100	-	ASTM D113